


VIRUS PRODUCTION PLATFORM BASED ON ADHERENT VERO CELLS

One platform for production of many viruses –
no repetitive media development.





VERO ADHERENT CELL LINE

The cell line is obtained from Public Health of England and adapted to serum free adherent growth. Fully adapted adherent cells are frozen in serum free, animal component free media. Cell bank history and development report are part of cell bank delivery package for high regulatory acceptance.



LOTUS – VERO ADHERENT CULTURE MEDIUM

LOTUS medium is free of serum, animal components and hydrolysates. Adherent VERO cells have been carried for more than 2 months in LOTUS medium with no loss of viability.

LOTUS medium contains recombinant proteins thus fulfill the highest regulatory standards.




Figure 2: Growth performance of VERO cells in commercially available media. Cells are adapted to corresponding media for 3 passages before testing the growth performance. Experiment is performed in T75 flasks and split ratio kept constant as 1:3.

Florabio AS
Gülbahçe Mah. Gülbahçe Cad.
Teknopark İzmir No: 01/43/04
35430 Urla-İzmir
TURKEY

Phone +90 232 502 3321
info@florabio.com.tr
www.florabio.com.tr